


IVTAGs...SIMPLE & EFFECTIVE


IMPROVE **COMMUNICATION**
 IMPROVE **PATIENT CARE**
 IMPROVE **OUTCOMES**

Easy snap on action. Plastic tag/labels for your devices.

Midlines



PICCs / CVCs



To Use:

1. Remove TAG from storage container.
2. Use ball point pen or fine permanent marker to write on TAG designated areas....**let ink dry. Prep-resistant markers provided with each order.**
3. Wrap extended tab around desired tubing and press snaps together, securing TAG to tubing.
4. TAGs are usable for all types of devices specific to their design.

Key communication features:

- Date
- Indication of device type
- Dwell Days
- Inappropriate uses
- Device specific considerations

Staff Improvements:

- Communication
- Nursing Compliance
- Patient Safety

Port Labeling Features:

- Date inserted
- Inappropriate Solutions

Product Features:

- Durable Lasting Plastic
- 2 Sided Communication
- Low Cost Item
- No glue
- Easy Application
- Never Expires
- Latex Free

www.IVTAGS.com



No Latex



Do not use if package is damaged